



CLA MEDICAL
HUMANISING HEALTHCARE

Registration NO:
IVDC9656524-189138

Medical Device
AUTHORITY
MALAYSIA




RAPID BIOTEC™ **Adenovirus / RSV / Influenza A+B Combo Rapid Test Device** (Nasopharyngeal Swab)

INTENDED USE

The Adenovirus, RSV and Influenza A+B Combo Rapid Test Device is a rapid chromatographic immunoassay for the qualitative detection of Adenovirus antigen, Respiratory Syncytial Virus antigen and Influenza A and B antigens in nasopharyngeal swab specimens. It is intended to aid in the rapid differential diagnosis of Adenovirus, Respiratory Syncytial Virus viral, Influenza A and B infections.

Operation Diagram



Performance Characteristics

Analytes	Sensitivity	Specificity	Accuracy
Adenovirus Antigen	98.2%	98.2%	98.2%
Influenza A+B	92.7%	98.0%	95.6%
Flu A	99.0%	98.9%	98.9%
Flu B	97.7%	99.0%	98.6%

Materials Provided

- Test Devices
- Workstation
- Sterile Swabs
- Package Insert
- Specimen Collection Tubes with Extraction Buffer

Products Advantages

- Easy specimen collection
- Fast results in 15 minutes assay time
- Easy visual interpretation
- Simple operation
- High accuracy

Ordering Information

Cat.no.	Product	Specimen	Kit Size
D-ARID10	Adenovirus, RSV and Influenza A+B Combo Rapid Test Device	Nasopharyngeal Swab	10T Device

For professional in vitro diagnostic use only.

Distributor

CLA MEDICAL SDN BHD (1336459-V)

+6012 859 6459

www.clamedical.com.my
MDA Licence No.: MDA-5480-WDP124