



CLA MEDICAL
HUMANISING HEALTHCARE

Registration NO:
IVDC9823725-193620

Medical Device
AUTHORITY
MALAYSIA



RAPID BIOTEC™


Cardiac Troponin I Rapid Test Device

(Whole Blood/Serum/Plasma)

INTENDED USE

The Cardiac Troponin I Rapid Test Device (Whole Blood/Serum/Plasma) is a rapid chromatographic immunoassay for the qualitative detection of human cardiac Troponin I in whole blood, serum or plasma as an aid in the diagnosis of myocardial infarction (MI).

Operation Diagram



Performance Characteristics

Analytes	Sensitivity	Specificity	Accuracy
Cardiac Troponin I	97.6%	99.4%	99.1%

Materials Provided

- Test Devices
- Prefilled Buffer
- Safety Lancets
- Droppers
- Package Insert
- Alcohol Swabs

Ordering Information

Cat.no.	Product	Specimen	Pack Size
D-TROPD20	Cardiac Troponin I Rapid Test Device	WB/S/P	20T Device

For professional in vitro diagnostic use only.

www.clamedical.com.my

MDA Licence No.: MDA-5480-WDP124

Distributor

CLA MEDICAL SDN BHD (1336459-V)

+6012 859 6459