



CLA MEDICAL
HUMANISING HEALTHCARE

Registration NO:
IVDC 7209024-186896

Medical Device
AUTHORITY
MALAYSIA



RAPID BIOTEC™


Dengue NS1 Rapid Test Device

(Whole Blood/Serum/Plasma)


INTENDED USE

The Dengue NS1 Rapid Test Device (Whole Blood/Serum/Plasma) is a rapid chromatographic immunoassay for the qualitative detection of NS1 antigen of Dengue virus in human whole blood, serum or plasma as an aid in the diagnosis of Dengue infections.

Operation Diagram



When tests become positive



Performance Characteristics

Analytes	Sensitivity	Specificity	Accuracy
Dengue NS1	95.8%	96.1%	96.0%

Materials Provided

- Test Devices
- Prefilled Buffer
- Safety Lancets
- Droppers
- Package Inser
- Alcohol Swabs

Ordering Information

Cat.no.	Product	Specimen	Kit Size
D-DAGD20	Dengue NS1 Rapid Test Cassette	WB/S/P	20T

For professional in vitro diagnostic use only.

Products Advantages

- Easy specimen collection
- Fast results in 10 minutes
- Easy visual interpretation
- Simple operation
- High accuracy

Distributor

CLA MEDICAL SDN BHD (1336459-V)

+6012 859 6459

www.clamedical.com.my

MDA Licence No.: MDA-5480-WDP124