



CLA MEDICAL
HUMANISING HEALTHCARE

Registration NO:
IVDC9656524-189138

Medical Device
AUTHORITY
MALAYSIA




Syphilis Rapid Test Device (Whole Blood/Serum/Plasma)

INTENDED USE

The Syphilis Rapid Test Device (Whole Blood/Serum/Plasma) is a rapid chromatographic immunoassay for the qualitative detection of antibodies (IgG and IgM) to *Treponema Pallidum* (TP) in whole blood, serum or plasma to aid in the diagnosis of Syphilis.

Operation Diagram



Performance Characteristics

Analytes	Sensitivity	Specificity	Accuracy
Syphilis Rapid Test Device	99.9%	99.7%	99.8%

Materials Provided

- Test Devices
- Safety lancets
- Package Insert
- Droppers
- Prefilled buffer
- Alcohol Swabs

Products Advantages

- Easy specimen collection
- Fast results in 5 minutes
- Easy visual interpretation
- Simple operation
- High accuracy

Ordering Information

Cat.no.	Product	Specimen	Pack
D-TPD20	TP (Syphilis) Rapid Test Device	WB / S / P	20 T

For professional in vitro diagnostic use only.